

# Bi-Elevation™

## sMAP3D™ FITTING GUIDE

Bi-Elevation™ is an advanced scleral alignment design feature used to help reduce instances of scleral misalignment


### 1 VERIFY THE LENS GEOMETRY AND DIAMETER

#### 1A. Determine lens diameter

- If HVID  $\neq$  11.7 mm, choose the 14.8-mm or 16.0-mm diameter lens
- If HVID  $>$  11.7 mm, choose the 15.4-mm or 17.0-mm diameter lens

#### 1B. Choose lens shape

- Choose **prolate** for eyes with central elevations
- Choose **oblate** for eyes with peripheral elevations or central flattening




1A

1B

## 2

## DETERMINE THE AMOUNT OF BI-ELEVATION™ TO APPLY TO THE LENS


**2A. Select "tool"**  to display toricity

**2B. Set "Chord Length"** to the initial landing point of your desired lens diameter

- The toricity at this chord length will determine the amount of "Bi-Elevation™" to add to your lens SAG (in this case 273 microns of Bi-Elevation™)

2A

2B



14.8 = 12.2 mm

16.0 = 12.8 mm

15.4 = 12.8 mm

17.0 = 13.6 mm

## 3


## SELECT YOUR FINAL SAG AND APS VALUES

**3A. Ensure lens diameter is correct** from step 1A

**3B. Set Shallow SAG based on SAG** in "Lens Parameters" (in this case 4.535mm)

- Add the Bi-Elevation™ amount from step 2B to this value for the second meridian SAG

**3C. Subtract the Bi-Elevation™ amount from the total toricity** (3C bubble) to determine the toricity in the APS



**Questions? Our expert consultants can help.** Gain individualized support from our expert fitting consultants available Monday to Friday, 8:00 AM to 7:00 PM EST. Call (800) 253-3669 or email [svp.consultation@Bausch.com](mailto:svp.consultation@Bausch.com)

Visit [bauschsvp.com](https://bauschsvp.com) for Important Safety Information.

®/™ are trademark of Bausch & Lomb Incorporated or its affiliates.  
Any other products/brand names and/or logos are trademarks of the respective owners.  
©2023 Bausch & Lomb Incorporated or its affiliates. ALZN.0032.USA.23

**BAUSCH + LOMB**  
**Zenlens**  
scleral lenses